

NIPT Sample Collection Instructions


Please refer to **Generation NIPT Collection SOP** for detailed instructions

Check the Form

- 1 Check that the gestational age is recorded and is more than 10 weeks at collection date
- 2 Check that the patient has read the consent information and has signed the request form
- 3 Check that the doctor has signed the request form

Collect and Mix Sample

- 1 Draw 7-10mL of blood into the supplied **Streck tube** using the Vacutainer system
A minimum of 7mL MUST be collected
- 2 **Immediately mix the tube by gently inverting 10 times**



- 3 Label the tube with the patient's full name and DOB in **DD/MM/YYYY** format
- 4 Record and initial the blood draw date in the Collection Information section of the form

Pack the Sample

- 1 Place blood tube, absorbant sheet and request form inside a Biohazard bag
- 2 Place the Biohazard bag inside the Genomic Diagnostics addressed self-seal envelope
- 3 Keep at room temperature until courier collection - **Do not refrigerate or freeze**

Generation NIPT Infoline: 1800 822 999

Generation
a new era in prenatal testing

Specialist Diagnostic Services Pty Ltd (ABN 84 007 190 043) t/a Abbott Pathology PUB/MR/1340_V6_Nov20

 Abbott Pathology